

On the instructions of a retained client

Engineering, veterinary, medical and laboratory equipment, etc

(Subject to availability)



Please be aware that this is a PPE site – visitors must wear adequate PPE at all times

ONLINE AUCTION:	To bid please go to www.lsh.co.uk/mba or  BidSpotter.co.uk <small>powered by atgmedia</small>
AT:	Building 819, Unit E, Discovery Park, Sandwich, CT13 9NJ
ON VIEW:	Monday 16 th August 2021, 10am – 4pm – <i>by appointment only</i>
BIDDING ENDS:	Tuesday 17 th August 2021 at 3pm
CLEARANCE:	Monday 23 rd – Friday 27 th August 2021, 9am – 4.00pm – <i>by appointment only</i>
FURTHER INFO:	Adam Marx amarx@lsh.co.uk
BID ENQUIRIES:	MBA Accounts 023 8046 1643 mbaaccounts@lsh.co.uk
OFFICE ADDRESS:	DeskLodge House, Redcliffe Way, Bristol, BS1 6NL

NOTICE TO PURCHASERS

VIEWING:

- The lots may be viewed as follows: **STRICTLY BY APPOINTMENT ONLY – Monday 16th August 2021, 10am to 4pm**

Please e-mail your viewing appointment request to amarx@lsh.co.uk.

- Located as follows: Building 819, Unit E, Discovery Park, Sandwich CT13 9NJ

BIDDING ENDS:

- Bidding will end at: **3pm on Tuesday 17th August 2021** [subject to the normal 10 min extension]
AT WHICH TIME ALL HIGHEST BIDS WILL BE CONSIDERED.

ACCEPTANCE OF FINAL HIGHEST BIDS – IMPORTANT NOTE:

- Acceptance of all final highest bids is subject to approval by our client.
- Lambert Smith Hampton reserves the right to reject any bids which they feel does not reflect reasonable value
- Successful bidders will be notified by email by **10am on Wednesday 18th August 2021.**

INSTRUCTIONS FOR BIDDING:

- Access for bidding can be obtained through our website www.lsh.co.uk/assetadvisory/online-auctions
- Full details including valid email address must be provided in order to receive approval for bidding
- In the event of a bid being received within 10 minutes of the scheduled closing time, the bidding period on that lot will automatically extend by a further 10 minutes and for a further 10mins for any bid received thereafter.
- **MAXIMUM BIDS:** If more than one maximum bid is placed of the same amount, the maximum bid placed first will be considered the highest bid and a further bid will need to be placed to take the lead. Bidders are notified by an automated message which appears on their screen
- **IMPORTANT:** When bidding, do not forget to 'REFRESH' your screen to ensure you are aware of current bids.
- **MAILING:** All approved bidders will automatically have their details placed on our database for notification of futures sales undertaken by Lambert Smith Hampton only. In the event that you wish to have your details excluded from our database, please use the 'unsubscribe' feature within our website (as above) alternatively please notify our Sales Administration Team aaaccounts@lsh.co.uk

BUYERS PREMIUM:

- The purchaser will pay a buyers premium of **16%** plus VAT, which will be charged on all lots.
- The premium is not negotiable and will be paid by all buyers

PAYMENT TERMS:

- All purchases are to be paid in full, no later than **2pm on Wednesday 25th August 2021**
- All queries regarding payment should be directed to our Sales Accounts team, aaaccounts@lsh.co.uk
- Please note that we do not have the facility to accept debit or credit card payments (the details required when signing up are for name and address validation purposes only).
- No cash will be accepted on sites for this sale – bank transfers only.

Please refer to the page headed "PAYMENT DETAILS" for acceptable methods of payment

PURCHASERS WHO MAKE OVER THE COUNTER PAYMENTS AT THE BANK WILL HAVE THEIR COLLECTION STALLED UNTIL THE BANK HAVE CONFIRMED IT IS A LEGITIMATE CASH PAYMENT.

VALUE ADDED TAX:

- Purchasers will be charged Value Added Tax on all lots, where applicable, at the current rates

DEPOSITS:

- Purchasers are required to pay a deposit on lots which will cause damage or leave a hole in the fabric of the building after removal. These lots will be specified in the sale catalogue. The deposit will be returned once the purchaser has made good in a satisfactory manner.

RISK:

- The bidder / buyer is at risk once the buyer is notified of his successful purchases and is strongly advised to effect insurance at once, irrespective of whether title has passed. Title does not pass to the purchaser until payment has been received in full.

RECORDED "USEAGE" READINGS: The hour and odometer readings on machinery and vehicles in this sale cannot be warranted and should not be relied upon.

COMPUTER SOFTWARE: Due to the Data Protection Act, all software and data will be removed from computers prior to selling.

CLEARANCE TERMS AND ARRANGEMENTS

- All lots must be removed between **STRICTLY BY APPOINTMENT ONLY – Monday 23rd – Friday 27th August 2021**
- The premises will be open from **9.00am to 4.00pm. No further access will be available after this time.**
- **All successful purchasers are to contact the auctioneer to arrange timeslots for collection. Numbers onsite during this period will be restricted due to the COVID-19 epidemic. Anyone who arrives without a prior collection appointment will be refused access to the site.**
- **Please note: There are no lifting/loading facilities on site.**
- **Please note: Scot Lift Systems Ltd can supply you with a free no obligation quotation to pack and ship your items won anywhere in the world. Please contact John W Walmsley on 01304 268010, 07788 443610 or johnwalmsley@scotlift.co.uk**

for more details. (Please note that payment of lots will need to be made directly to Lambert Smith Hampton prior to delivery. Lambert Smith Hampton does not guarantee the services provided by any third party service provider and accepts no responsibility for any problems experienced. All arrangements are between yourselves and the third party service provider)

- Lambert Smith Hampton does not undertake packaging, postage, delivery or shipment of goods and the bidder must make his own arrangements to comply with the clearance terms
- Clearance of all lots must be undertaken in accordance with Health & Safety at Work Regulations and, where necessary, Construction Design and Management Regulations 1994.
- **Risk Assessments and Method Statements (RAMS):** If you purchase a lot that requires RAMS to be completed please be aware that you will not be permitted access to site to commence collection of ANY of your lots until the required RAMS has been completed and reviewed by an LSH representative. If you need access to site during the collection period to review equipment for RAMS purposes, or have any further questions please contact Adam Marx, amarx@lsh.co.uk.
- Lambert Smith Hampton reserve the right to stop clearance if in their opinion, they believe it is being carried out in an unsafe manner or without the Method Statement and Risk Assessment.
- The purchaser should ensure that their contractors hold sufficient public liability insurance, copies of which should be made available to Lambert Smith Hampton if demanded.
- Any fluids remaining in any lots purchased MUST be removed from site in conformity with the Control of Substances Hazardous to Health (COSHH) Regulations.

CONDITIONS OF SALE

All lots will be sold subject to the Conditions of Sale as printed on the catalogue cover. Bidders are strongly advised to read them carefully.

TRANSFER OF LOTS:

Lambert Smith Hampton will only accept payment from and permit removal of lots purchased by the successful purchaser or their duly appointed agent.

OVERSEAS PURCHASERS:

Overseas purchasers must, on acceptance of offers, advise of their intentions to export their good immediately and notify our Sales Co-ordinator, Elaine Allan [eallan@lsh.co.uk]

Overseas purchasers will be required to provide details of :

- (a) Their passport(s)
- (b) One other method confirming their identity, ie drivers licence, name, address, UK agent (if any), telephone and fax numbers.
- (c) Full details of their national bank and their UK corresponding bank.
- (d) Details from their bank confirming method of payment, etc.

Please also ensure that:

- (e) Your country, or the country to which the items are destined, holds no import restrictions on the goods that you wish to purchase.
- (f) The country has no import licence restrictions or a restriction on currency allocation.
- (g) The country has no prolonged inspection procedure, which might cause excessive delay in allowing your goods to be imported.
- (h) You receive a fully descriptive invoice in order that you may arrange payment as soon as possible with our bank.
- (i) You employ a reputable freight-forwarding organisation. Please note that many of the machinery removal companies in this country are not necessarily freight forwarders. We will be pleased to advise accordingly.
- (j) You have read the Conditions of Sale.

Please Note:

- (a) We will not allow removal of any Lot from the Sale site until our Accounts Department has received notification from our bankers that your payment has been properly cleared in full. Monies must be transferred to ourselves within the time stated in the "Notice to Purchasers"
- (b) If there are any special arrangements that you wish to make with us, you must contact us at least three days before the Sale date.
- (c) Our terms regarding payment of VAT and subsequent refunds.

VAT on Goods Destined for Export:

All purchasers intending to export their purchases will be charged a deposit equal to that of the VAT due on each lot purchased and the Buyer's Premium.

In order to qualify for a refund, the goods must be exported in accordance with regulations under "Indirect Export" and valid evidence of export received by Lambert Smith Hampton within one month of the date of export. Where satisfactory evidence is not provided, goods and services will be subject to UK VAT at the standard rate and the deposit will not be refunded.

Please ensure that the cashier has been notified of your intentions to export prior to removal of goods.

The above information must be supplied in writing to our Sales Co-ordinator, Mrs Elaine Allan, Lambert Smith Hampton, First Floor, 2 Manor Court, Barnes Wallis Road, Fareham, PO15 5TH, on your official company notepaper and signed by the contact name quoted.



If you intend to export your purchases using your own transport (either by sea, rail or air), then VAT will be levied and refunded upon receipt of satisfactory documentation, supporting removal of the goods purchased.



Lambert Smith Hampton reserve the right to refuse to accept the bid of any bidder should they be unable to satisfy Lambert Smith Hampton on any of the above.

Scotlift Systems Limited (LabAssets)

Address:

Building 819, Unit E, Discovery Park, Sandwich, CT13 9NJ





NOTICE

**ALL LOTS IN THIS
SALE ARE SUBJECT
TO A BUYERS
PREMIUM OF 16%
PLUS V.A.T AT THE
PREVAILING RATE**

PAYMENT DETAILS

◆ Purchasers should note that **GOODS MAY ONLY BE REMOVED FROM SITE UPON RECEIPT OF PAYMENT IN FULL IN CLEARED FUNDS**

◆ For the avoidance of doubt, CLEARED FUNDS is:

Telegraphic Transfer / BACS [see note 2]

Credit / Debit cards / Cash / Bankers Drafts / Cheques CANNOT be accepted

PURCHASERS WHO MAKE OVER THE COUNTER PAYMENTS AT THE BANK WILL HAVE THEIR COLLECTION STALLED UNTIL THE BANK HAVE CONFIRMED IT IS A LEGITIMATE CASH PAYMENT.

(2) Telegraphic Transfers / BACS payments	Purchasers who wish to make direct payments into our bank account can obtain bank details from the office following notification of successful bids. Please note timescales for payment, as follows <ul style="list-style-type: none">◆ Telegraphic payments : up to 24hrs◆ BACS payments : 3 day Payments direct into our bank account will not be accepted by any other method. For the avoidance of doubt over the counter credits are not an acceptable method of payment and is in breach of Lambert Smith Hampton payment terms particularly those that relate to the Money Laundering Act 2017

MONEY LAUNDERING REGULATIONS

Services provided by Lambert Smith Hampton in relation to the sale of assets fall under the Proceeds of Crime Act 2002. Accordingly we have put in place policies and procedures which are designed to prevent, forestall and identify acts of money laundering.

Lambert Smith Hampton will not accept transactions in cash which exceed £7,000 under any circumstances and Lambert Smith Hampton will take any action necessary to ensure that it does not contravene its obligations under the Money Laundering Regulations 2017. In the event that a purchaser entering into a business relationship with Lambert Smith Hampton, knowingly or unwittingly fails to comply with our policy Lambert Smith Hampton will have no choice than to treat this as a breach of the regulations. Where necessary ID will be sought by the company and consider a report to NCIS (National Criminal Intelligence Service)

Should you have any queries on the above, please contact our AA Accounts Team on +44(0)2380 461643

Lot No.	Description
---------	-------------

1	Vacuum Drying Oven, DZF-6125 complete with operating manual
2	Grant JB Aqua 2 max 37.5 Deg.
3	Graphtec Thermal Arraycorder WR8500 (WA13279)
4	Thermo Scientific S1 Pipette Filler with charger
5	GE Power Protection Unit M-Pact Plus Air Circuit Breaker up to 4000a - unused
6	Metal Racking system which slots together 34 uprights 1980mm x 600 127 x 1200mm Shelf supports 57 x 900mm Shelf supports 33 x 1800mm Shelf supports 13 x 1500mm Shelf supports 126 Shelves 58mm x 59mm 41 Shelves 89mm x 47mm
7	27 Bays of warehouse racking, 24 bays blue and 3 Bays Black Each Bay 1000mm x 490mm deep x 2140mm High 4 shelves in each bay plus a roof panel allowing you to store.
8	E2M18 Oil Sealed Rotary Vacuum Pump.(WA13710)
9	Two PE Sciex Line adjustment option (WA10953)
10	Harvard Apparatus Medium Hydraulic Surgical Table, SWL 150kgs.(WA13430)
11	Gould 5900 Signal Conditioner Cage Model CL-810231-01, design to facilitate production of a family of products and provides signal conditioner outputs via individual BNC connectors or Bare Wire connectors. The 5900 has been designed to meet and tested to international safety and EMI/RFI standards. (WA11248_
12	WTW Multi 340i Set. (WA13014)
13	Oxford Cryosystems controller 600 series (WA11097)
14	Industrial Vacuum system 3 Phase
15	APEX mill Serial no. M 2070 (Ex Pfizer so well maintained)
16	Tank system can be used with Lot 15 (Ex Pfizer so well maintained)
17	Argo Ocean twin 1 Aircon unit and condenser
18	Cedex Innovatis GmbH Automated Cell Counting System. (WA12867)
19	Carl Zeiss HTS-Reader 340 nm. (WA12868)
20	Panasonic KX-B430 write on electronic board
21	SympaTec GmbH Rodos and Quixel are instruments used in analysing wet/dry, particle/droplet size and shape. The RODOS is an universal dry dispersing unit applicable for nearly all kinds of dry powders with particle sizes in the range from < 0.1µm to about 3500 µm. Quixel is the advanced wet dispersing system suitable for particle size analysis of all kinds of suspensions and

**Engineering, veterinary, medical and laboratory equipment
Bidding Ends: Tuesday 17th August 2021 at 3pm**

Lot No.	Description
---------	-------------

	emulsions in the size range of 0.1 µm up to 3.5 mm. (WA13200)
22	Mettler Toledo platform balance - visually compromised
23	Caportec Protein Detector
24	UVD 340U is a high performance photodiode array detector for use with modular HPLC systems. The design of the detector is optimized for maximum sensitivity in the UV range, while maintaining sufficient resolution to distinguish spectra even of closely related compounds. (WA11253)
25	Various Thermo Scientific Pipettes and stand (WA10987)
26	FitzMill Pump - Not Checked (WA11083)
27	Alcatel Pascal 2015 SD Dual Stage Rotary Vane Vacuum Pump. (WA13784)
28	Sarel 83318 Electric cabinet, IP66, 400mm x 300mm x 200mm
29	Two 25ml Hex base glass measuring Cylinders FB55205
30	Shoreline underbench Sparkfree Refrigerator
31	GE Power Protection Unit M-Pact Plus Air Circuit Breaker up to 4000a - Unused
32	Sarel 83318 Electric cabinet, IP66, 400mm x 300mm x 200mm
33	BVC PP188 Process Vacuum System (WA13342)
34	Mettler Toledo Automatic filling machine
35	Sieve 300mm
36	Foramaflow VAP/2 Formalin Vaporiser, 2000ml. (WA13352)
37	Avery Berkel KL24 Model L117/L217/L227 Intrinsically Safe Weighing System
38	Avery Berkel KL22 Model L117/L217/L227 Intrinsically Safe Weighing System
39	Tall Freezer with Locking door (WA11114)
40	Three Spear and Jackson Garden Rakes - unused
41	CopyDisc label printer and robotic arm (WA11784)
42	Lab Products Enviro-guard Environmental Control System Supply Unit, Model 59013. (WA12819)
43	Lab Products Enviro-guard Environmental Control System Supply Unit, Model 59016N (WA12818)
44	Water Aspirator SH2-95B Circulating Water Vacuum Pump New in Box can be used with large Evaporator for extraction - unused
45	Biotage Isolera LS system Biotage® Isolera LS is the definitive flash chromatography system for reliable scaling up from grams to multi, even hundred gram scale purification, using the industry standard original Isolera interface. Enjoy all of the sophisticated software features of laboratory scale systems in a robust instrument built for scale-up labs. Isolera Spektra LS adds an integrated pump assisted loading mechanism for larger samples, and a funnel rack collection option which allows samples to be collected in larger fraction collection vessels.

**Engineering, veterinary, medical and laboratory equipment
Bidding Ends: Tuesday 17th August 2021 at 3pm**

Lot No. Description

46	The iCELLis nano lid forms a part of the iCELLis bioreactor system, which is an automated, single-use, fixed-bed bioreactor that provides cell growth conditions for adherent cells.
47	Manual Furniture Movers 500kgs loading <i>NB: This item has no record of Thorough Examination. The purchaser must ensure a Thorough Examination is carried out by a competent person prior to first use. The purchaser will be required to sign an undertaking confirming their agreement</i>
48	Roll of steel banding still in original Packaging
49	Taylor Wharton HC35 dewar is used for storing samples in liquid nitrogen as well as storing and dispensing small amounts of liquid nitrogen. Features include:- Static holding time - 130 days- Working time - 81 days- Normal Evaporation rate (litres/day) - 0.27- Liquid Nitrogen Capacity - 35 litres- Vial Capacity 1020 (WA13654)
50	Sarel 83318 Electric cabinet, IP66, 400mm x 300mm x 200mm
51	Penton RCS 6/T Recessed Metal Ceiling Speaker - unused in box
52	One-channel digital indicator displays measured values and monitor limit values. The universal input allows for the connection of thermocouples, RTD temperature probes, resistance transmitters, potentiometers, or measuring probes with standard signal. The detection rate for the measuring input is 4 measurements per second. The display indicates the measured values and allows dialog control during configuration. (WA13230)
53	Leica Friomobile Freezer Unit. (WA12817)
54	Sartorius TE612 Talent Top-loading Analytical Balance, 610 g x 0.01 g (WA11903)
55	iCELLis nano lid forms a part of the iCELLis bioreactor system, which is an automated, single-use, fixed-bed bioreactor that provides cell growth conditions for adherent cells. (WA13408)
56	Thermo Electron Corporation orion 3 star ph benchtop single-channel meter with calibration and measurement accuracy. Reliable, accurate pH measurement for any general laboratory or QC application. (WA11836)
57	Graduated and one-mark Easypet hand-held pipette that responds exactly to finger pressure for slow or rapid aspiration, setting the meniscus or emptying the pipette: Suitable for pipettes of 0.1-100 mL. High aspirating speed: 25 mL in 5 sec. Cordless work for up to seven hours, continuing directly after recharging. Autoclavable pipette adapters and filter holders. Active pressure compensation prevents corrosion caused by aspirated vapour. unused in Box (WA12950)
58	Large Gantry lifting Beam with trolley 3200mm x 850mm x 2230mm High unbolts for transport <i>NB: This item has no record of Thorough Examination. The purchaser must ensure a Thorough Examination is carried out by a competent person prior to first use. The purchaser will be required to sign an undertaking confirming their agreement</i>
59	Wheaton W348923-D bench top mini bottle roller apparatus for cell cultivation research or staining procedures. Caters for serum

**Engineering, veterinary, medical and laboratory equipment
Bidding Ends: Tuesday 17th August 2021 at 3pm**

Lot No.	Description
---------	-------------

	and media bottles, as well as culture tubes. (WA11846)
60	Masterflex Precision Pump Tubing 25FT
61	Thermomax C0428 SM Quattro (-50 / +50) microprocessor-based datalogger. (WA12807)
62	Wooden Export packing case with latches, approx. size 900mm x 800mm x 1020mm high
63	Packing Case 1380mm x 880mm x 1400mm high
64	BD FACSCelesta flow cytometer is designed to make multicolor flow cytometry more accessible and allow researchers to benefit from new innovations in instrument and reagent technology. Parallel advances in optical and reagent technology are pairing multilaser, multidetector instruments with bright new dyes, enabling increasingly deep and powerful insights in cell analysis, complete with all hardware and software plus many additional items. (WA13576)
65	BOC Medical Regulator, HP1800S, with 2x Stewart Buchanan Pressure Gauges, max outlet 800psig. (WA13342)
66	Clinitek Status+ Analyser from Siemens is an easy-to-use, automated analyser designed to deliver accurate and reliable results using Siemens urinalysis strips, provides a comprehensive point-of-care urinalysis menu for routine testing, detection of early kidney disease and hCG pregnancy testing. Unused (WA13686)
67	Fisherbrand Digital Colorimeter Model 45 unit measures the absorbance of particular wavelengths of light by a specific solution and filters. General Specification:- Resolution, 0.02 at 1Abs (absorbance), 2% (transmission).- Wavelength, 440 to 700nm.- Range, 0 to 199% (transmission).- Operating Temperature Range, 0 to 199% (transmission).- Bandwidth, 40nm.
68	Generic Glass Test Tubes, 12x75mm, single box. (WA13447)
69	Ashcroft 'use no oil' 0-600 psi gauge on Premier Industries 50-10916-01 regulator, Inlet: 500 PSIG/Flow: 0.5 SLPM, and proprietary brass low pressure gas regulator, with Rodent Face Masks and Feeder Tubes, non-rebreathing system type.
70	Penton RCS 6/T Recessed Metal Ceiling Speaker - unused in box
71	Long Read Tower automated DNA sequencer based on VGI's MicroGene Clipper platform. (WA10629)
72	Spirax Monnier IR1M 3.5 bar compressed air regulator for general purpose pneumatic systems. (WA13341)
73	Large Range of Stainless Steel Surgical Instruments
74	PALL Life Systems The iCELLis nano lid forms a part of the iCELLis bioreactor system, which is an automated, single-use, fixed-bed bioreactor that provides cell growth conditions for adherent cells. (WA13415)
75	Sarel 83318 Electric cabinet, IP66, 400mm x 300mm x 200mm

**Engineering, veterinary, medical and laboratory equipment
Bidding Ends: Tuesday 17th August 2021 at 3pm**

Lot No.	Description
---------	-------------

76	Thermo Scientific REVCO REL 404V solid door high performance laboratory refrigerator is ideal for storing samples that have clinical, industrial and scientific applications. It uses a microprocessor control system, positive airflow systems, industrial-quality cabinet construction and extra strength refrigeration compressors.(WA11697)
77	U S Gauge Model 23 15psi, max inlet 400psig, and Contents 400 psi gauge gas regulators with rodent face masks, non-rebreathing system type. (WA13308)
78	Penlon MRI IM 500 CO2/N2O anaesthetic gas regulator machine. (WA11711)
79	Alpha unit block assembly for PTC DNA engine systems (WA10806)
80	Multidrop 384 Reagent Dispenser continuous high-speed dispensing of various liquids with precision. It can portion a wide range of solutions into 96- and 384-well plates across a 5 to 395L volume range, providing reliable microvolume dispensing for pharmaceutical and biotechnology laboratories. (WA13590)
81	Sarel 83318 Electric cabinet, IP66, 400mm x 300mm x 200mm
82	Fedegari Autoclave FVA/A1 vertical lab sterilizer, Model FVA2/A1, is equipped with a stainless steel chamber and lid with automatic independent operation and control. The steriliser has been specially designed with ergonomical features for easy handling of sterilization treating materials and goods including: glassware, liquids in open and sealed containers, mixed wastes/regular and high risk, as well as surgery instruments, textiles, rubber articles for hospital use. The model has the counter pressure sterilization option (steam-air mixture for LVP, liquids in sealed containers and heat sensitive products).
83	Silogic Veterinary EC60 Portable Cardiac and Respiratory Monitor. (WA11764)
84	One-channel digital indicator displays measured values and monitor limit values. The universal input allows for the connection of thermocouples, RTD temperature probes, resistance transmitters, potentiometers, or measuring probes with standard signal. The detection rate for the measuring input is 4 measurements per second. The display indicates the measured values and allows dialog control during configuration.
85	Elga Inline Water Filter System (WA11718)
86	Biodex 069-310 surface monitor is compact 3-range surface rate meter, equipped with an internal 2-inch diameter pancake GM detector and built-in speaker. The unit is an appropriate tool for surveying bench tops and other surfaces, and for checking hands and clothing for radioactive contamination. The meter face reads in both cpm and mR/h. Anti-saturation circuitry keeps meter needle at full scale in high radiation fields. Gamma Sensitivity: Nominal is 3000 cpm/mR/h (based on Cs 137). (WA13709)
87	Rigaku X-Stream XS-2000 Crystallography System, with Cold Head and Nozzle and including Control Box and Helium Compressor. The X-stream 2000 Cryogenic Crystal Cooler is a self-contained apparatus designed for cooling crystals to a

**Engineering, veterinary, medical and laboratory equipment
Bidding Ends: Tuesday 17th August 2021 at 3pm**

Lot No.	Description
---------	-------------

	temperature within the range of -180C to +4C.
88	Sorvall Heraeus 4 swing bucket rotor body with four place swing bucket rotor, no buckets. Max 5850rpm. (WA13092)
89	DiLab AccuSampler Standard automated Blood Sampling (ABS) (WA11176)
90	Columbus Instruments Micro Injector 500 (WA11088)
91	Sarel 83318 Electric cabinet, IP66, 400mm x 300mm x 200mm
92	Multiskan Spectrum Verification Plate, Catalog No. 240 73405, and case Catalog No. N03394.
93	Penlon Loud Speaker, unused
94	DiLab AccuSampler is a fully automated system for blood sampling (ABS) in freely moving laboratory rodents such as mice, rat, guinea pig, hamster. The system handles sample volumes between 25 and 500 ul. (WA12695)
95	Assorted Glass bundle including many different items of Glass with a minimum of 40 Items. Could contain evaporator glassware
96	Kontron Micromon 7142 Vital Signs Monitor with Microrecorder 7341 Chart Recorder. (WA13264)
97	Penton RCS 6/T Recessed Metal Ceiling Speaker - unused in box
98	Integra Biosciences Vacusafe Aspiration system is a compact, all-in-one laboratory vacuum pump, providing a safe and efficient way to collect and contain biological liquid waste. (WA13054)
99	Eppendorf F45-30-11 30-place fixed-angle basic rotor with PTFE-coated body and lid. Capacity: 30x 1.5/2.0ml tubes, Max speed: 20,817 x g (14,000 rpm). (WA12945)
100	Nova Biomedical System fluid cartridge Bioprofile Flex, 39580. (WA13398)
101	Teledyne Isco RediSep Alumina Neutral functionalized columns are designed to help to improve column consistency and reliability. (WA12994)
102	Cytomics FC 500 Flow Cytometer Automated tube-based acquisition for clinical assays. The system has the capacity to conduct 5-color analysis using either a single (488nm) or a dual (488nm/635nm) laser configurations. (WA13584)
103	VWR flat wire gauze, with ceramic centre, provides excellent support to beakers and other glassware during heating on ring stands. (WA13652)
104	10mm bladed hole borer in hand made wooden box
105	Two 25ml Hex base glass measuring Cylinders FB55205
106	One-channel digital indicator displays measured values and monitor limit values. The universal input allows for the connection of thermocouples, RTD temperature probes, resistance transmitters, potentiometers, or measuring probes with standard signal. The detection rate for the measuring input is 4 measurements per second. The display indicates the measured

**Engineering, veterinary, medical and laboratory equipment
Bidding Ends: Tuesday 17th August 2021 at 3pm**

Lot No.	Description
---------	-------------

	values and allows dialog control during configuration.
107	Amersham Pharmacia Biotech LeadSeeker Homogeneous Imaging System with 2 x Filter Control Unit. (WA13713)
108	Fisherbrand Digital Colorimeter Model 45 unit measures the absorbance of particular wavelengths of light by a specific solution and filters. General Specification:- Resolution, 0.02 at 1Abs (absorbance), 2% (transmission).- Wavelength, 440 to 700nm.- Range (WA12172)
109	Eppendorf F45-30-11 is a 30-place fixed-angle basic rotor with PTFE-coated body and lid. Capacity: 30x 1.5/2.0ml tubes, Max speed: 20,817 x g (14,000 rpm). (WA12944)
110	Seven SKF Roller Bearings 6006-2RS1
111	Fisher Scientific FB61222 Total Range Thermometer Traceable (WA13201)
112	Digital Colorimeter Model 45 unit measures the absorbance of particular wavelengths of light by a specific solution and filters. General Specification:- Resolution, 0.02 at 1Abs (absorbance), 2% (transmission).- Wavelength, 440 to 700nm.- Range
113	Penton RCS 6/T Recessed Metal Ceiling Speaker - unused in box
114	Mettler-Toledo 51340288 InLab 413 SG Combination pH/ATC Electrode. (WA11755)
115	Teledyne Isco RediSep Alumina Neutral functionalized columns are designed to help to improve column consistency and reliability. (WA12993)
116	Dionex Universal Chromatography Interface UCI-100 Universal Interface, for use with ion chromatography. Chromatographic detectors convert signals. Data collection rate (100 Hz) and high channel stability. Four independent RS-232 interfaces. (WA11261)
117	BioHit Proline 8 Channel Electronic Pipette 50-1200ul. (WA13001)
118	BioHit Proline 8 Channel Electronic Pipette 50-1200ul. (WA13002)
119	Thermo Electron Corporation orion 3 star ph benchtop, single-channel meter with improved stability algorithm to increase calibration and measurement accuracy, accurate pH measurement for any general laboratory or QC application.
120	U S Gauge Model 23 15psi/Contents 400psi Gas Regulators with Rodent Face Masks and Feeder Tubes (WA13308)
121	Vetronic Services Ventilator SAV03 Small Animal Ventilator Controller. (WA13268)
122	Teledyne Isco RediSep Alumina Neutral functionalized column, designed to help improve column consistency and reliability. (WA12991)
123	Buddeberg Laboratory Technology PLR 11 Compressed Air Laboratory Stirrer, 4026446002510, designed for a variety of mixing tasks. The vane type, air-driven motor offers a high power

**Engineering, veterinary, medical and laboratory equipment
Bidding Ends: Tuesday 17th August 2021 at 3pm**

Lot No.	Description
---------	-------------

	density and is significantly lighter and smaller than an equivalent electric motor. (WA13316)
124	Pall Emflon II, 5 inch Cartridge Gas Filter, 0.2um, SLK7002V002PV. (WA13403)
125	BioHit Proline 8 Channel Electronic Pipette 25-250 uL (WA12457)
126	Consort EV243 Electrophoresis Power Supply. General features:- Timer 0-99.59 hour.- Volt hours 0-99.99 kVH.- LCD display, 2 16 characters.- Ground leakage detection 500 uA.- No load detection, programmable, overload detection.- 3600 values of data logging.- 0-95% non-condensing humidity.- 3 4 mm sockets, 30 ohms min load resistance.- limit = 0-40 C temperature.- output = 0-400 V, 0-300 mA, 0-50 W. (WA11811)
127	Teledyne Isco RediSep Alumina Neutral functionalized columns are designed to help to improve column consistency and reliability. (WA12992)
128	Generic Glass Test Tubes, 12x75mm, single box. (WA13447/2)
129	Hawksley Haemacytometer Set
130	Penton RCS 6/T Recessed Metal Ceiling Speaker - unused in box
131	Consort EV231 Electrophoresis Power Supply. General features:- Timer 0-99.59 hour.- Volt hours 0-99.99 kVH.- LCD display, 2 16 characters.- Ground leakage detection 500 uA.- No load detection, programmable, overload detection.- 3600 values of data logging.- 0-95% non-condensing humidity.- 4 4 mm sockets, 70 ohms min load resistance.- limit = 0-40 C temperature.- output = 0-1,200 V, 0-500 mA, 0-300 W. (WA11808)
132	Fisherbrand Digital Colorimeter Model 45 unit measures the absorbance of particular wavelengths of light by a specific solution and filters. General Specification:- Resolution, 0.02 at 1Abs (absorbance), 2% (transmission).- Wavelength, 440 to 700nm.- Range, 0 to 199% (transmission).- Operating Temperature Range, 0 to 199% (transmission).- Bandwidth, 40nm. (WA12172)
133	Consort EV231 Electrophoresis Power Supply. General features:- Timer 0-99.59 hour.- Volt hours 0-99.99 kVH.- LCD display, 2 16 characters.- Ground leakage detection 500 uA.- No load detection, programmable, overload detection.- 3600 values of data logging.- 0-95% non-condensing humidity.- 4 4 mm sockets, 70 ohms min load resistance.- limit = 0-40 C temperature.- output = 0-1,200 V, 0-500 mA, 0-300 W. (WA11807)
134	BioHit Proline Electronic Pipette Charging Carousel. (WA12499)
135	Dispensette Digital Easy Calibration Organic Bottle top dispenser
136	Biohit Proline 8 Channel Electronic Pipette 250. (WA12458)
137	Sartorius Biohit 791202 1200?L Non-sterile Polypropylene Optifit Tips, Refills. Pack of 960. (WA13675)

**Engineering, veterinary, medical and laboratory equipment
Bidding Ends: Tuesday 17th August 2021 at 3pm**

Lot No.	Description
---------	-------------

138	Jencons Sealpette Electronic Pipette Charging Carousel. (WA12500)
139	One-channel digital indicator displays measured values and monitor limit values. The universal input allows for the connection of thermocouples, RTD temperature probes, resistance transmitters, potentiometers, or measuring probes with standard signal. The detection rate for the measuring input is 4 measurements per second. The display indicates the measured values and allows dialog control during configuration.
140	Eight boxes of Thermo Scientific Frosted Microscope Slides, Ground 90. (WA13657)
141	Two 25ml Hex base glass measuring Cylinders FB55205
142	Gilson 306 claims reliable chromatographic separation, utilizing easy method development through TRILUTIONLC software or using the stand-alone keypad. The Interchangeable pump heads accommodate flow rates ranging from 5 L/min - 200 mL/min with a wide range of liquids at differing viscosities. Pumps are available in isocratic, binary, or tertiary configurations. (WA13565)
143	Dionex Universal Chromatography Interface UCI-100 Universal Interface, for use with ion chromatography. Chromatographic detectors convert signals. Data collection rate (100 Hz) and high channel stability. Four independent RS-232 interfaces. (WA11258)
144	Alpha Unit single-block assembly, model 384, holds one 384-well microplate or one 96-well microplate. (WA12890)
145	Gast Pump. Model 1531/107b/g289x (WA10813)
146	Multiposition Microelectric Valve Actuator Type EMHMA-CE. Four valves. (WA11500)
147	Sarel 83318 Electric cabinet, IP66, 400mm x 300mm x 200mm
148	Pfizer bench top waste liquid pump (WA13504)
149	Bio-Rad GS-800 calibrated densitometer, system 95-800 for gel, film and blot analysis. It has a built-in calibrated optical transmissive step tablet for accuracy. Each tablet contains a serial number for completing documentation to an external standard, ensuring accurate and reproducible results. (WA11678)
150	Buddeberg PLR 11 compressed air laboratory stirrer is designed for a variety of mixing tasks. The vane type, air-driven motor offers a high power density and is significantly lighter and smaller than an equivalent electric motor. (WA13319)
151	Jasco PTC-351S Peltier Controller. (WA13041)
152	Gas-Arc Lab-Master Hydrogen, CO ₂ , Oxygen & Nitrogen Pressure Regulators (13 off) and CC Isolation Va (WA12924)
153	Bioseb Rodent Pincher RP1 Algometer. (WA13043)
154	Jasco DG_1580_32 dynamic mixer is a high pressure dynamic mixer for mixing up solvents under high pressure. (WA11283)

**Engineering, veterinary, medical and laboratory equipment
Bidding Ends: Tuesday 17th August 2021 at 3pm**

Lot No.	Description
---------	-------------

155	MSD Sector Imager 6000 demonstration plate, 036ZCAI1188. (WA12938)
156	HEL Heating Block (WA11521)
157	Wesbart IS89 Multi Shelf Shaker Incubator. (WA13040)
158	Masterflex Precision PumpTubing 25FT
159	Whirlimixer (WA11472)
160	Model 21 30psi, max inlet 500psig, gas regulators. (WA13301)
161	Alpha unit block assembly for PTC DNA engine systems (WA10809)
162	Hawksley Haemocytometer Set
163	IKA Color Squid small magnetic stirrer. Electronically controlled motor for more capacity. High speed range. Max. stirring quantity 1 l. Outstanding chemical resistance with glass top and synthetic base. (WA11724)
164	Stainless steel storage box with stainless steel lid and tray, and with polished surface finish. (WA13461)
165	Waters Pump Control Module II PC2 (WA11726)
166	Vetronic Services ERM-8010 Veterinary ECG and Respiratory Monitor. (WA13269)
167	Two 25ml Hex base glass measuring Cylinders FB55205
168	Honey Comb Fraction Collector (WA10976)
169	Two Multiposition Microelectric Valve Actuator Type EMTMA-CE (WA11501)
170	Replacement charcoal absorber for closed cycle gas compressors. (WA13489)
171	Test tubes are manufactured from borosilicate glass with rimless wall design, ideal for chemical and biological work. (WA13256)
172	Type 832 Multidrop 384 dispensers is a well microplate liquid dispensing system which is compact, easy to use and maintain. Features that prevent cross contamination and save on reagents. For higher throughput needs, Multidrop dispensers can be easily integrated to robotic setups. (WA12448)
173	Fluostar Optima and Polarstar Optima Service Box, BSAN-500/G-752. (WA13302)
174	Sarel 83318 Electric cabinet, IP66, 400mm x 300mm x 200mm
175	Analytical Sales Temperature Control Module, model PTC050, Hot Sleeve Column Heater is designed to offer optimum performance in heating applications. It has fast 'H-Bridge' circuitry that claims to improve control responsiveness. Its variable time base burst firing increases the heater life and gives better temperature controllability. (WA13272)
176	DeGasys DG-2410. Two channels degasser uses fluororesin membrane for gas transmission and separation. Flowrate 7 ml/min. per channel. Liquid contact parts: PTFE, PEEK, PPS . Dead volume per channel: 10 ml. (WA11729)

**Engineering, veterinary, medical and laboratory equipment
Bidding Ends: Tuesday 17th August 2021 at 3pm**

Lot No.	Description
---------	-------------

177	KQ Integrated Solutions APCI Probe OPTON-20012, P/N 97055-60090 Rev A - unused in box (WA13273)
178	Fisherbrand Digital Colorimeter Model 45 unit measures the absorbance of particular wavelengths of light by a specific solution and filters. General Specification:- Resolution, 0.02 at 1Abs (absorbance), 2% (transmission).- Wavelength, 440 to 700nm.- Range, 0 to 199% (transmission).- Operating Temperature Range, 0 to 199% (transmission).- Bandwidth, 40nm. (WA12174)
179	Philips Focusline Projection, White Light (90-100 Ra), Type 5761 30W 6V Halogen Lamp. (WA12931)
180	GE Healthcare Whatman HEPA-Cap 75 Disposable Air Filtration Capsule, 6702-7500. (WA13394)
181	GelAir Drying System, for drying polyacrylamide and agarose gels. Use in densitometry, photo documentation, autoradiography, overheads, and long-term storage. The heated drying chamber works like a convection gel ovens though drying times may vary depending on percentage and thickness of the gel. (WA12842)
182	Pall Emflon II, 5 inch Cartridge Gas Filter, 0.2um, SLK7002V002PV. (WA13396)
183	Penlon Sigma Delta Key Fill Isoflurane Vaporizer - unused in box (WA12834)
184	Hawksley Haemacytometer Set
185	Generic Glass Test Tubes, 12x75mm, single box. (WA13447/3)
186	Assorted Glass bundle including many different items of Glass with a minimum of 40 Items. Could contain evaporator glassware
187	Econoneb is a compressor that is ideal for both the hospital and the clinic environments. It claims a high flow rate that is suitable for nebulising bronchodilators (relievers), corticosteroids (preventers) and antibiotics. (WA13507)
188	IKA Color Squid small magnetic stirrer. Electronically controlled motor for more capacity. High speed range. Max. stirring quantity 1 l. Outstanding chemical resistance with glass top and synthetic base. (WA11723)
189	Five DrySyn Wax bowls 50ml (WA13160)
190	Mitutoyo Digital Micrometer 0-1inch, 293-832. (WA13179)
191	Two 25ml Hex base glass measuring Cylinders FB55205
192	Sartorius Sartopore 2 XLG Sterile Capsule. Pore size: 0.8 0.2 m (WA13393)
193	Dispensette Digital Easy Calibration Organic Bottle top dispenser
194	Drager Pac III claims to provide a fast and accurate response to hazardous gas levels using an audible alarm to alert the user that the danger is present. The chrome-plated housing also provides resistance to potential Radio Frequency Interference (RFI). (WA13508)

**Engineering, veterinary, medical and laboratory equipment
Bidding Ends: Tuesday 17th August 2021 at 3pm**

Lot No.	Description
---------	-------------

195	Dorset ID Identification Trovan LID1260 DESKTOP READER, version 1260V724m, is a desk reader that can read all Trovan transponder types with a LCD 2 x 16 characters with backlight, viewing area of 99 x 24 mm. EEprom for data storage, 5400 slots without data/time stamp, 2900 with data/time stamp. (WA13500)
196	OmniPAGE VS10DCAST 10 x 10cm casting base unit, precast gel unit for routine mini protein electrophoresis using Mini 10x10cm vertical gel systems. (WA12219)
197	Critikon Dinamap 1846 SX Vital Signs Monitor with Integrated Oxytrak Pulse Oximeter (WA13254)
198	One-channel digital indicator displays measured values and monitor limit values. The universal input allows for the connection of thermocouples, RTD temperature probes, resistance transmitters, potentiometers, or measuring probes with standard signal. The detection rate for the measuring input is 4 measurements per second. The display indicates the measured values and allows dialog control during configuration.
199	Grass Instrument Regulated Power Supply, Model RPS 107E. (WA13255)
200	Hawksley Haemocytometer Set
201	Penton RCS 6/T Recessed Metal Ceiling Speaker, unused in box
202	Reaction Glass Vessel - excluding framework
203	Alpha unit block assembly for PTC DNA engine systems (WA10808)
204	Thermos 4 litre Container, Type 3504. (WA12555)
205	Jasco DG-2080-53(3) on-line HPLC vacuum degasser removes gaseous components such as O ₂ , N ₂ , and CO ₂ from a solvent by passing the solvent through a fluoropolymer membrane tube with the pressure outside the tube reduced. Vacuum degassing by membrane separation is carried out under mild conditions and has no effect on the composition of the solvent. General specification: - Maximum flowrate: 10ml/min (with water). - Maximum usable flowrate: 4ml/min. - Maximum pressure: 2.0kgf/cm ² . - Solvent hold-up volume: Approx. 10ml each line. (WA11282)
206	REMS Autosampler (WA11079)
207	Anelva Thermocouple Gauge TG-550 - unused in box
208	Roche lightcycler S/N 1402740 Used for PCR Testing (WA11190)
209	Assorted Glass bundle including many different items of Glass with a minimum of 40 Items. Could contain evaporator glassware
210	Agilent 1100 Series WP G2255A Robotic Well-plate Handler Stacker. (WA12947)
211	ETTAN DALT II Gel Caster, Model 80-6467-22. (WA13048)
212	Hawksley Haemocytometer Set
213	Genetix AliQuot ALQ Reagent Dispenser Microplate Plate Washer. (WA13186)

**Engineering, veterinary, medical and laboratory equipment
Bidding Ends: Tuesday 17th August 2021 at 3pm**

Lot No.	Description
---------	-------------

214	Alpha unit block assembly for PTC DNA engine systems (WA10807)
215	Heidolph MR 3004 safety C Magnetic Stirring Hotplate (WA13326)
216	Dilvac Medium and Small Wax Flasks, 2pcs. (WA13088)
217	Philips Focusline Projection, White Light (90-100 Ra), Type 5761 30W 6V Halogen Lamp. (WA12931)
218	Two 25ml Hex base glass measuring Cylinders FB55205
219	Reaction Glass Vessel - excluding framework
220	Hawksley Haemacytometer Set
221	Four DrySyn Wax bowls 100ml (WA13163)
222	Ecoline VC-380 V2.13 - Microprocessor controlled tubing pump (WA11323)
223	VistaLab Technologies Ovation Bionatural Macro Adjustable Volume Electronic Pipette, 02ul to 10ml. (WA12705)
224	Prior Precision Microscope Stage Type HK01E600. (WA11629)
225	One-channel digital indicator displays measured values and monitor limit values. The universal input allows for the connection of thermocouples, RTD temperature probes, resistance transmitters, potentiometers, or measuring probes with standard signal. The detection rate for the measuring input is 4 measurements per second. The display indicates the measured values and allows dialog control during configuration.
226	PULSAFE Respiratory Protection. (WA11238)
227	Binocular microscope - 10x eyepieces with integral objective lens: 6/ 12/ 25/ 50/. (WA11590)
228	DIONEX UVD 340U High Performance Photodiode Array Detector. (WA11253)
229	Millipore water purifier Milli Q plus (WA11243)
230	Reaction Glass Vessel
231	IPGphor 80-6414-02 IEF Isoelectric Focusing Unit performs first dimension isoelectric focusing (IEF) in 2-dimensional (2-D) electrophoresis. The IPGphor features programmable functions including rehydration time, platform temperature, maximum current limit, voltage change pattern, voltage for each step in the isoelectric focusing program, and step duration. Unit comes with case of IPGphor strip holders (6pcs) (WA13378)
232	Remp manual piercer (WA11008)
233	Rheodyne Labpro Valve Controller, model PR100-106-02. (WA11001)
234	DIONEX Universal Chromatography Interface UCI-100. (WA11262)
235	Five DrySyn Wax bowls 100ml (WA12822)
236	Alpha unit block assembly for PTC DNA engine systems (WA10806)

**Engineering, veterinary, medical and laboratory equipment
Bidding Ends: Tuesday 17th August 2021 at 3pm**

Lot No.	Description
---------	-------------

237	Reaction Glass Vessel - excluding framework
238	Amersham Pharmacia Biotech IPGphor Strip Holder Electrophoresis (2pcs). (WA13379)
239	Mettler-Toledo 51340288 InLab 413 SG Combination pH/ATC Electrode. (WA11755)
240	Hawksley Haemocytometer Set
241	IKA Dispensing Tool 100-2000ML, S 25 N - 25 F - unused (WA13226)
242	Jenway 3320 pH Meter. (WA11032)
243	Vitalograph BREATHCO Carbon Monoxide Monitor Model 29700, 2pcs. (WA13299)
244	Variomag Thermomodule 40ST Controller with Magnetic Stirring Hotplate. (WA13185)
245	Jasco DG-2080-53(3) on-line HPLC vacuum degasser removes gaseous components such as O ₂ , N ₂ , and CO ₂ from a solvent by passing the solvent through a fluoropolymer membrane tube with the pressure outside the tube reduced. Vacuum degassing by membrane separation is carried out under mild conditions and has no effect on the composition of the solvent. General specification: - Maximum flowrate: 10ml/min (with water).- Maximum usable flowrate: 4ml/min.- Maximum pressure: 2.0kgf/cm ² .- Solvent hold-up volume: Approx. 10ml each line. (WA11281)
246	Locking Brackets for Oxford AS 500 NMR (WA13365)
247	Jasco CO-1560 intelligent column thermostat (WA11287)
248	Instron Impact Accessories. (WA12922)
249	Rheodyne Labpro Valve Controller, model PR100-106-02. (WA11000)
250	Colmark C95031S Pressure meter (WA10842)
251	IKA Color Squid (WA11014)
252	Leica 600 mm Cold Light Source CLS Goose Neck, 30130011 - unused (WA13208)
253	DeGasys DG-2410. Two channels degasser uses fluororesin membrane for gas transmission and separation. Flowrate 7 ml/min. per channel. Liquid contact parts: PTFE, PEEK, PPS . Dead volume per channel: 10 ml. (WA11277)
254	Gast Pumps, Series 1531 (WA10853)
255	JENCONS whirlmix vx 100 S/N;71221056 (WA11199)
256	Carl Zeiss light source (WA10965)
257	IKA Color Squid (WA11015)
258	Novex PowerEase 500 is designed specifically for mini-gel electrophoresis. It offers extensive programming capabilities, including eight pre-set and four custom methods for you to set to your own preferences. The simple, intuitive PowerEase interface is faster and easier to use than traditional dial-in power supplies. In addition, PowerEase features: constant voltage, current, or

**Engineering, veterinary, medical and laboratory equipment
Bidding Ends: Tuesday 17th August 2021 at 3pm**

Lot No.	Description
---------	-------------

	power settings; buffer temperature monitoring capability; on the fly program editing; and compatibility with most printers using a parallel port. (WA11038)
259	EMKA Technologies Ventilator Vent2. (WA13220)
260	Stem Corporation RS1000 Block Heater. Compact ten position (well) heating and stirring reaction block. Temperature up to 150C. Stir speed 400-2000 rpm. External control capability via RS32/RS4851 interface. Suitable for integration onto a robotic platform. (WA11002)
261	Avisoft Bioacoustics UltraSoundGate 116 (WA13087)
262	Solid Stainless Steel and Glass Reaction Vessel.
263	Multidrop Micro Accessories. (WA13056)
264	One-channel digital indicator displays measured values and monitor limit values. The universal input allows for the connection of thermocouples, RTD temperature probes, resistance transmitters, potentiometers, or measuring probes with standard signal. The detection rate for the measuring input is 4 measurements per second. The display indicates the measured values and allows dialog control during configuration.
265	Aeroqual R10 Remote Sensor Head Adapter. The Aeroqual R10 remote sensor head adapter lets you measure gases in one area, while you remain safely in another. It is designed for applications where the sensor needs to be separate from the monitor base. (WA13053)
266	No lot
267	Eiko AV/Photo lamp, EFR 15v 150w is a halogen MR16 lamp with a GZ6.35 two-pin base. It features a rated colour temperature of 3400K and an average lifetime of 50 hours. 2pcs. (WA13058)
268	GE EJA 21V 150W Replacement Halogen lamp for Projectors and Microscopes, 2pcs. (WA13468)
269	Gilson 306 chromatographic separator, utilizing easy method development through TRILUTIONLCsoftware or using the stand-alone keypad. The Interchangeable pump heads accommodate flow rates ranging from 5 L/min - 200 mL/min with a wide range of liquids at differing viscosities. Pumps are available in isocratic, binary, or tertiary configurations.
270	Drummond Portable Pipette-Aid XL, with a longer than normal handle, performs the same operation as many other conventional pipettors. The ergonomic design includes an adjustable sliding hand rest. The Portable Pipette-Aid XL also features a quieter, more powerful pump.(WA13140)
271	In Vitro Systems & Services Vivascience Universal Turning Device. (WA13569)
272	BMDS Data Acquisition System, DAS-6006 Probe and IPTT 4365 Link Model DAS-7007S. (WA13070)
273	Stuart Scientific Auto Vortex Mixer SA2. (WA12420)
274	Assorted Glass bundle including many different items of Glass

**Engineering, veterinary, medical and laboratory equipment
Bidding Ends: Tuesday 17th August 2021 at 3pm**

Lot No.	Description
---------	-------------

	with a minimum of 40 Items. Could contain evaporator glassware
275	Assorted Centrifuge Accessories. (WA13116)
276	EMD Millipore, Merck Molecular sieve 0.3 nm, 1kg, 1057041000. (WA13653)
277	Matrix Technologies Impact Electronic Multichannel Pipette, 850ul, with stand and power (WA13120)
278	Radleys Discovery Technologies Greenhouse Parallel Synthesiser & Blowdown Evaporator (WA13114)
279	LINK fourth order (four-pole) noise filter to improve detection in liquid chromatography. It has been developed as a stand-alone unit and is applicable with any detector that supplies an analogue signal. The improvement in Signal-to-Noise ratio depends on the frequency of the signal vs. the noise. (WA12410)
280	Biomat 2 biological cabinet with stand
281	Stainless Steel and Trespa bench with castors
282	WIKA Duel Gauge Medical Gas Regulator. (WA13068)
283	Hawksley Haemocytometer Set
284	Hanna pHep 1 Waterproof Pocket-sized pH Tester with 0.1 pH Resolution H198151 (WA13069)
285	DeGasys DG-2410. Two channels degasser uses fluororesin membrane for gas transmission and separation. Flowrate 7 ml/min. per channel. Liquid contact parts: PTFE, PEEK, PPS. Dead volume per channel: 10 ml. (WA11275)
286	Burkard Scientific Benchtop Safety Cabinet. (WA13190)
287	Millipore Stainless Steel Process Valve (WA10801)
288	Assorted Glass bundle including many different items of Glass with a minimum of 40 Items. Could contain evaporator glassware
289	BMDS Data Acquisition System, DAS-6006 Probe and IPTT 4365 Link Model DAS-7007S. (WA13071)
290	One-channel digital indicator displays measured values and monitor limit values. The universal input allows for the connection of thermocouples, RTD temperature probes, resistance transmitters, potentiometers, or measuring probes with standard signal. The detection rate for the measuring input is 4 measurements per second. The display indicates the measured values and allows dialog control during configuration.
291	UGO Basile 5025 Dog Ventilator, 700 ml cylinder/ piston. It consists of a high-precision motor operated by a solid state speed control, coupled to a 700 ml cylinder/piston assembly via an original variable stroke linkage. The mechanism does not incorporate any slotted link. Tidal volume and frequency can be adjusted while the ventilator is running or standing still. The piston almost touches the cylinder end with each stroke, regardless of the preset volume. This insures that all air taken into the pump is expelled with each stroke.

**Engineering, veterinary, medical and laboratory equipment
Bidding Ends: Tuesday 17th August 2021 at 3pm**

Lot No.	Description
---------	-------------

292	Edwards CG16K capsule dial gauge is barometrically compensated with NW flange fittings. Designed to cover a range of 0 to 1040 mbar, the gauge provides accurate, repeatable performance even at low pressures making them equally suited for non-corrosive process plant or for laboratory applications. (WA13381)
293	Fume Cube Purex Soldering Fume Extraction System
294	Motic binocular microscope body - 10x/21 eyepieces with integral objective lens: 0/6/ 12/ 25/ 50/. (WA12527)
295	Two 25ml Hex base glass measuring Cylinders FB55205
296	Stainless Steel and Trespa bench with castors
297	IPGphor 80-6414-02 IEF Isoelectric Focusing Unit performs first dimension isoelectric focusing (IEF) in 2-dimensional (2-D) electrophoresis. The IPGphor features programmable functions including rehydration time, platform temperature, maximum current limit, voltage change pattern, voltage for each step in the isoelectric focusing program, and step duration. Unit comes with case of IPGphor strip holders (6pcs) (WA13377)
298	OmniPAGE VS10DCAST 10 x 10cm casting base unit, precast gel unit for routine mini protein electrophoresis using Mini 10x10cm vertical gel systems. (WA12221)
299	Two 25ml Hex base glass measuring Cylinders FB55205
300	Marble Balance slab
301	Stainless Steel and Trespa bench with castors
302	Trespa Laboratory Bench with wheels, resistant to all chemicals and can withstand well over a tonne of loading (WA13743)
303	DentiMex MiniTech EL 750 Dental Hand piece Drill Unit. (WA13573)
304	GE Power Protection Unit M-Pact Plus Air Circuit Breaker up to 4000a - unused

**Engineering, veterinary, medical and laboratory equipment
Bidding Ends: Tuesday 17th August 2021 at 3pm**

1: Conditions

These conditions, together with any 'Special Conditions' or 'Notices to Purchasers', are the only terms and conditions subject to which Lambert Smith Hampton (LSH) (the "Auctioneer") acting as agents for the Vendor (the "Vendor") will contract with buyers (the "Buyer"). Bidding for any lot shall be deemed to be an acceptance of these conditions and any conditions in the Notices to Purchasers. If there is any inconsistency between any of these conditions and anything in the Notices to Purchasers, then those in the Notices to Purchasers shall prevail.

2: Inspection and Description

- 2.1 Buyers have a responsibility to make their own inspection and investigation of the lots at the time(s) and at the premises specified, and to satisfy themselves on all matters affecting the lots, and to inspect and satisfy themselves prior to the sale, as to the condition and description of a lot, its fitness and suitability for purpose.
- 2.2 Goods are believed to be correctly described but all goods are sold with faults, imperfections and errors of description.
- 2.3 Neither the Vendor nor the Auctioneers, their servants or agents makes or gives, nor has any person in the employment of the Auctioneers any authority to make or give, any representation or warranty in relation to any lots.

3: Limitation of Liability

- 3.1 Neither the Vendor nor the Auctioneers, their servants or agents shall be liable for any loss or damage suffered by the Buyer arising out of or in connection with any defects or deficiencies in any lots purchased, errors of description in the auction catalogue, any mis-statements as to any matter affecting the lots or the failure of the lots to fulfill the functions for which they were intended
- 3.2 The Auctioneer's and the Vendor's total liability for breach of contract or negligence is limited to the Bid Price.
- 3.3 Neither the Auctioneers nor the Vendor shall be liable for any indirect or consequential loss or damage (whether for loss of profit, loss of business or otherwise) incurred by the Buyer whatsoever.
- 3.4 Nothing in these conditions excludes or limits the liability of the Auctioneers or Vendor for death or personal injury caused by their negligence or fraudulent misrepresentation.
- 3.5 Nothing in these conditions shall prejudice the statutory rights of a consumer who shall be bound by these conditions only insofar as they are consistent with such statutory rights.

4: Conduct of Sale

- 4.1 The Auctioneers reserve the right in their absolute discretion to refuse admission to the premises and grounds to any person.
- 4.2 The conduct of the sale shall be at the Auctioneer's sole discretion and the Auctioneers reserve the right to refuse any bids without giving any reason. If any dispute of whatsoever nature relating to a bid shall arise, such dispute shall be determined by the Auctioneers in their absolute discretion and whose decision shall be final and binding on all parties.
- 4.3 Each bidder must register prior to submitting any bids in Online Sales. The Auctioneer reserves the right to reject any registration at its sole discretion.
- 4.4 The Vendor and the Auctioneers reserve the right to fix a reserve price for any lot.
- 4.5 No bid may be withdrawn after the close of bidding.
- 4.6 The Vendor or the Auctioneers have absolute discretion to withdraw consolidate or divide at any time any lot.
- 4.7 The highest bidder for each lot at the close of bidding shall be the Buyer, subject to approval by the Auctioneers and the Vendor. The Sale Price is exclusive of the buyer's premium and VAT. The highest bid at the close of bidding subject to approval by the Auctioneers and Vendor also marks a conclusion of a contract of sale between the Vendor and the Buyer.
- 4.8 In making a bid for any lots, the Buyer does so as principal. The Buyer will be held personally and solely liable for a bid unless the Auctioneers have previously agreed in writing with the Buyer that the bidder will do so on behalf of an identified third party acceptable to the Auctioneers. In circumstances where the Auctioneers have so agreed, the Buyer and the third party will jointly and severally be liable for all obligations arising from the bid where a third party shall be bound by these conditions by the individual bidding as his agent in the same way as if he were bidding personally.

5: Price

The price payable by the Buyer for the purchase of the lots secured shall be the sale price, being the highest bid and accepted by the Auctioneer in accordance with condition plus the buyers premium and any other associated costs (if any) and VAT applicable at the prevailing rate.

6: VAT

The price of the goods purchased is exclusive of Value Added Tax (VAT), which where applicable, the Buyer shall be additionally liable to pay to the Auctioneer

7: Risk and Title

- 7.1 Insurance: As from the date and time the Buyer is declared, all risks in and relating to the Goods purchased, shall pass to the Buyer and the Buyer is strongly advised to effect insurance for these risks at once. In no circumstances will the Auctioneers or the Vendor be responsible if any Goods or part thereof be lost, stolen, damaged or destroyed after the Buyer is declared.
- 7.2 Title: The Vendor shall only sell such title to the Goods as he may have. Neither the Auctioneers nor the Vendor warrants the Vendor's good title to any of the Goods and if it is found that the Vendor does not have title or unencumbered title to any of the Goods purported to be sold under these conditions the Buyer expressly agrees that it shall have no right either to rescind the contract or to claim damages or a reduction in the consideration paid or payable under the contract.
- 7.3 Title to each of the Goods will not pass until:-
 - (a) all debts owed by the Buyer to the Auctioneers (whether part-paid secured or otherwise) are settled, for the avoidance of doubt the debts shall not be settled until (where payment or part payment is made by cheque) all cheque(s) in question have been cleared; and
 - (b) such Goods have been removed from the premises at which the sale is held or where the Goods are being stored in their entirety.

8: Removal

- 8.1 No lot or part thereof can be removed until payment in cleared funds (cash, approved Bankers Draft, BACS or Chaps) has been made in full or (where any part of the payment was made by cheque) the cheque has cleared in the Auctioneer's account. No lot shall be removed without the authority of the Auctioneer and unless under the supervision of the Auctioneer's servants or agents.
- 8.2 Removal shall take place only up to 5.00 pm and must be completed no later than the dates and times specified and published within the Notice to Purchasers or any date and time specifically declared thereafter, in writing, by the Auctioneer or written authority is obtained from the Auctioneers for clearance after that time.
- 8.3 The Buyer has no right to anything not described in the auction catalogue.
- 8.4 Delivery of lots sold will be made only to the declared Buyer and no transfer of any lot or part of a lot to any other person will be recognised.
- 8.5 Prior to the removal of any lot the Vendor may rescind the contract for the sale of that lot and refund to the Buyer any money paid by the Buyer for the lot should any third party claim title to or possession of any part of the lot.
- 8.6 The Buyer will be responsible for removal of lots at his own expense and must provide his own labour, and equipment.
- 8.7 The Auctioneers require that in pursuit of safe working practice all equipment used for lifting and transportation should be covered by appropriate insurance, safety documentation set out by current Health & Safety legislation at the time and registration documents (in particular, but not limited to, cranes, lifting equipment, tackle and fork lift trucks). Such documentation shall be produced to the Auctioneers and if the Buyer fails to produce such documentation on request for inspection the Auctioneers reserve the right to refuse to release the lot and/or rescind the sale in which event the provisions of condition 8 will apply. Equipment may only be used by operatives who have the valid and appropriate licences to use such equipment
- 8.8 Appropriate Personal Protective Equipment (PPE) must be worn during these activities and due regard paid to safe material handling practices.

9: Damage to Premises

The Buyer shall be responsible for the removal of all goods purchased at his own expense and risk and such removal must be carried out safely and lawfully and in accordance with any conditions of sale of which the Buyer is notified. For the avoidance of doubt removal includes where applicable disconnection from the mains electricity supply and the detachment of any lots which are fixed. All electrical installations must be carried out by a Certified electrician and left in a safe condition during and after the removal and the use of explosives or flame cutting equipment or any other potentially hazardous or inflammatory process shall not be permissible on the site without the express written consent of the Auctioneers. The Buyer shall be responsible to insure against and to make good any injury or damage to persons or property caused by the Buyer, their carriers, servants or agents detaching, disconnecting or removing any goods purchased by the Buyer. The Buyer shall produce his insurance documentation to the Auctioneers on request and if the Buyer fails to produce such documentation the Auctioneers reserve the right to refuse to clear the lot and/or rescind the sale. The Buyer shall indemnify the Auctioneers for any damage or loss which the Auctioneers may suffer in respect of loss, damage or injury suffered by the Buyer's agents or employees or any third party arising from the removal. At its discretion the Auctioneers may require the payment of a deposit to the Auctioneers prior to and/or during the removal which will be refunded on the removal being completed to the Auctioneer's satisfaction. If the Buyer refuses to deposit such monies, the Vendor or the Auctioneers may refuse the Buyer access to the premises for the purposes of removing the lot purchased by them and/or rescind the sale in which event the provisions of condition 8 will apply. The Auctioneers reserve the right to stop the removal of a lot by a Buyer, his agents or employees, if it appears to the Auctioneers or its agents or employees that the removal is being carried out in an unsafe or unsatisfactory way.

10: Default by Buyer

Upon failure by the Buyer for whatsoever reason:-

- (i) to pay a deposit in full if required under the Conditions of Sale or special notices; or
- (ii) to pay the purchase price in full on the due date; or
- (iii) to provide any documentation required under these conditions; or
- (iv) to remove the goods in a safe or satisfactory way; or
- (v) to remove any lot from the premises on or before the date specified for removal,

The Vendor will be entitled to rescind the contract forthwith without incurring any liabilities to the Buyer and, without prejudice to any claims of the Vendor and/or the Auctioneers against the Buyer arising from breach of contract or otherwise, upon rescission as aforesaid the following provisions apply:

- (a) all monies deposited in part payment will be forfeited and used to pay the Vendor's and/or Auctioneer's expenses referred to in paragraph (d) below;
- (b) if the lots have been removed in breach of these conditions, the Vendor and/or the Auctioneers, their servants or agents may enter the premises of the Buyer to recover such lots;
- (c) lots may be re-sold or otherwise disposed of by the Auctioneers in the manner they feel appropriate at their sole discretion and any deficiency arising upon the resale together with the expenses of it shall be due as a debt from the Buyer in default upon the first sale. The Buyer consents to such a resale on the Auctioneer's conditions of business applicable at the time of resale. Neither the Auctioneers nor the Vendor shall be liable to account to the Buyer in the event of a re-sale at a higher price than the price contracted to be paid by the Buyer. The Buyer waives any claim in such a case that he may have title to the lot and agreed that any re-sale price shall be deemed commercially reasonable;
- (d) the Buyer will be responsible for the specified time for final removal for all losses and expense incurred by the Vendor and/or the Auctioneers including storage, security and removal expenses, the costs of re-selling or disposing of lots and Auctioneer's commission.

11: Auctioneer as Agent

The Auctioneers act only as Agents for and on behalf of the Vendor and shall not be held responsible for any action or default on the part of either the Vendor, bidders or the Buyer. Any concluded contract of sale is made directly between the Vendor and the Buyer.

12: Accident or Damage

Neither the Auctioneers nor the Vendor will accept any responsibility for any accident, (except for that arising out of its negligence, or the

negligence of its agents or employees, resulting in death or personal injury) howsoever caused to any person which may occur whilst on the Auctioneer's premises, the Vendor's premises or such other premises used to hold the auction or to store the lots before during or after the sale and any person entering the premises does so at his own risk and is deemed to have notice of the condition of the premises and their contents.

13: The Consumer Protection Act 1987

No lots are sold as new.

14: Health & Safety At Work

It is expressly brought to the Buyer's attention and that of potential buyers, at the time of sale, any item of plant, machinery or equipment contained in the goods may not necessarily comply with the Health and Safety at Work Act 1974; the Environmental Protection Act 1990; Construction Design and Management Regulations 2015 or any other Act of Parliament or regulations thereunder governing the use of that plant, machinery or equipment in a working environment. Buyers of any such plant, machinery or equipment are hereby required to ensure so far as reasonably practicable that such item will be safe and without risk to health and that the use of any such items at a place of work within the EEC & the United Kingdom does not contravene any such Act of Parliament or regulation applicable to such use.

15: Dangerous Substances

It is expressly brought to the attention of Buyers (and potential buyers) that certain types of plant or main service installations could contain blue and white asbestos, dangerous chemicals and hazardous waste which if not handled correctly during their removal from the premises could be in breach of the Health and Safety at Work Act 1974, Sections 2-13 Control of Substances Hazardous to Health Regulations 2002 (COSHH) or other current legislation regulating the use of such substances in a working environment. The Buyer shall ensure so far as is reasonably practicable that they comply with the Health and Safety at Work Act 1974, COSHH and other current legislation in connection with the removal, handling and transport of such dangerous substances and or hazardous waste or shall employ a specialist contractor to remove them. The Buyer shall, on request, satisfy the Auctioneers in relation to its removal procedures, and the removal of waste materials must be undertaken by an approved and licensed contractor. The Buyer will indemnify the Vendor, their servants or agents against any loss, damages or expenses suffered by them as a result of the Buyer's or a sub-contractor's failure to comply with this legislation.

16: Notices

Any notices or other communications shall be in writing and, if sent by post, shall be deemed to have been received by the addressee on the second working day after posting or, if the addressee is outside the United Kingdom, on the fifth working day after posting. If any written notice is delivered by hand, it shall be treated as having been received at that time. Any notice sent to the Auctioneers shall be sent to the Auctioneer's address as set out in the auction catalogue, or displayed on the auctioneer's web site on the notice specific to that sale. Any notice which the Auctioneers send to the Buyer may be sent to the Buyer's last address known to the Auctioneer.

17: Third Party Rights

Save as expressly provided in these conditions, no term of these conditions shall be enforceable under The Contracts (Rights of Third Parties) Act 2001 by a third party.

18: Governing Law

These Conditions and any Conditions contained in the Notices to Purchasers, along with all associated transactions and all connected matters shall be governed and construed in accordance with English Law.

19: Severability

In the event that any provisions of these conditions shall be held unenforceable for any reason, the remaining conditions shall remain in full force and effect.